



Manual no. 106018

Item no.: Model

- 80131 VELA Latin 100
- 80132 VELA Latin 300
- 80111 VELA 'Basic' Sampling Chair



The chair's design varies depending on the model

TABLE OF CONTENTS

1.0. INTRODUCTION	3
1.1. SAFETY	4
1.1.1. Usage	4
1.2. WARRANTY	4
1.3. UNPACKING THE CHAIR	4
1.4. MAINTENANCE	5
1.4.1. Wheels	5
1.4.2. Frame	5
1.4.3. Upholstery	5
1.5. DISPOSAL	5
1.6. USER MANUAL	6
1.6.1. Attaching the backrest	6
1.6.2. Adjustment of seat height with foot ring*	6
1.6.3. Adjustment of seat height with lever	6
1.6.4. Adjustment of seat tilt	6
1.6.5. Adjustment of back height	7
1.6.6. Adjustment of back tilt	7
1.5.7. Footrest ring*	7
1.6.8. Armrest for sample collection*	7

1.0. INTRODUCTION

DEAR CUSTOMER

Congratulations on your new VELA Medical Chair. This user manual contains useful information about positioning, operating and maintenance.

Important

Please read this user manual thoroughly and keep it for later use. This user manual must always accompany the chair. The user manual is also available on VELA Medical's website: www.vela-medical.com, where you also can find other relevant information about the chair.

If you have any questions, please contact VELA Medical's distributor in your country.

VELA Medical reserves the right to make changes to this manual.

VELA Medical

Gøteborgvej 8-12
DK-9200 Aalborg SV
Denmark
www.vela-medical.com

GENERAL

1.1. SAFETY

This chair bears the CE mark and approved under the model name "VELA Latin" and conforms to all applicable EU rules. The chair is manufactured by VELA Medical, Goeteborgvej 8-12, 9200 Aalborg SV, Denmark.

1.1.1. USAGE

The chair is an indoor chair designed for sitting persons.

- :: The chair is not built to stand on.

The Chair is tested and approved for a maximum user weight of 125 kg.

1.2. WARRANTY

VELA Medical provides warranty by the laws and regulations which are debt in the country the chair is bought in. The warranty only applies when using original spare parts and accessories. Modification may only be carried out by authorised service technicians.

VELA Medical does not accept responsibility for damage or injury to the product or user, which have arisen due to:

- :: Transportation
- :: Misuse
- :: General wear and tear
- :: Incorrect use

- :: Use of non-original spare parts and accessories
- :: Adjustments carried out by unauthorised service technicians

1.3. UNPACKING THE CHAIR

The chair is supplied in cardboard packaging, developed with consideration for the environment.

If there are any parts missing or transport damage to the chair, please contact your dealer immediately. Do not try to repair any damage or use the chair.

Content of the cardboard packaging:

- :: VELA Medical chair
- :: This user manual

Note: see section. 1.6.1. for correct mounting of the backrest

1.4. MAINTENANCE

We recommend your chair undergoes an annual service to find, for example, defective or worn parts. The user should regularly retighten all removable parts.

1.4.1. WHEELS

VELA Medical recommend that the wheels are checked in regular intervals e.g. to remove possible thread remnants.

1.4.2. FRAME

Clean with water, possibly with a PH neutral soap.

1.4.3. UPHOLSTERY

Fabric: Clean with upholstery cleaners or wipe it with a tightly wrung cloth rinsed in clean water.

Imitation leather: Clean with water or mild soap.


Medical textiles: Clean as described above. Can also be cleaned with disinfectants and alcohol. To see which specific concentrations the different pad types can handle, see www.vela-medical.com

1.5. DISPOSAL

This product contains recyclable materials such as plastic, metal, etc. As a result, product disposal must be carried out in compliance with local legislation and separately from regular household waste. Correct disposal and recycling contribute to preventing the product from having a negative impact on the environment and people.

1.6. USER MANUAL

Replacement of the gas spring may only be performed by the service support team from your dealer or a skilled technical person.



1.6.1. ATTACHING THE BACKREST

The lever for backrest adjustment (1) must be in open position. Press the locking pin spring (2) and push the backrest arm through the hole. Tighten the lever afterwards.

Note: The backrest arm must be pushed through completely so the pin spring locks the backrest arm. Make sure that the backrest arm cannot be pulled out again.



1.6.2. ADJUSTMENT OF SEAT HEIGHT WITH FOOT RING*

Adjust the seat height with your foot by pushing the ring, which is placed at the base of the chair.



1.6.3. ADJUSTMENT OF SEAT HEIGHT WITH LEVER

Adjust the seat height by pushing the lever placed at your right hand.



1.6.4. ADJUSTMENT OF SEAT TILT

Adjust the seat tilt by pushing the middle lever placed at your right hand.



1.6.5. ADJUSTMENT OF BACK HEIGHT

Adjust the back height by using the lever placed on the back of the chair. Remember to tighten the handle again afterwards.



1.6.6. ADJUSTMENT OF BACK TILT

Adjust the back tilt by pushing the rear lever placed at your right hand.



1.5.7. FOOTREST RING*

A 45 cm footrest ring, adjustable in height. The footrest ring provides stability, support and comfort in the sitting position. A footrest ring is especially practical when the user is sitting high up.



1.6.8. ARMREST FOR SAMPLE COLLECTION*

The armrest is adjust by loosening the screw under the cushion and regulate the angle of the armrest. Remember to tighten the screw afterwards.

** Depending on model*

