



VELA Medical

Manual no. 106065

Item no.: Model

5027 VELA 'Bariatric' Patient Chair

5107 VELA 'Bariatric+' Patient Chair


ISO 9001
ISO 14001*The chair's design varies depending on the model*

TABLE OF CONTENTS

| | |
|--|-----------|
| 1.0. INTRODUCTION | 3 |
| 1.1. WARNINGS | 4 |
| 1.2. TROUBLESHOOTING | 5 |
| 1.3. VELA TANGO WITH ELECTRIC HEIGHT ADJUSTMENT | 6 |
| 1.3.1. Power input | 6 |
| 1.3.2. Emission | 6 |
| 1.4. LABELS USED ON THE CHAIR | 6 |
| 1.5. SAFETY | 7 |
| 1.5.1. Intended use | 7 |
| 1.5.2. Transport | 7 |
| 1.6. WARRANTY | 8 |
| 1.7. UNPACKING THE CHAIR | 8 |
| 1.8. MAINTENANCE | 8 |
| 1.8.1. Wheels | 8 |
| 1.8.2. Frame | 8 |
| 1.8.3. Padding | 8 |
| 1.8.4. Recycling guide | 9 |
| 1.9. RECYCLING AND DISPOSAL | 9 |
| 2.0. USER MANUAL | 10 |
| 2.0.1. Height adjustment (El)* | 10 |
| 2.0.2. Height adjustment (Manual)* | 10 |
| 2.0.3. Handheld control (El)* | 10 |
| 2.0.4. Adjustment of seat depth | 10 |
| 2.0.5. Adjusting the backrest angle | 11 |
| 2.0.6. Setting the back height | 11 |
| 2.0.7. Adjusting the armrest height | 11 |
| 2.0.8. Adjusting the armrest width | 11 |
| 2.0.9. Brake* | 12 |
| 2.0.10. Charging (El)* | 12 |
| 2.0.11. Strolling bracket* | 12 |
| 2.0.12. Electric brake (El)* | 12 |
| 2.0.13. 4-wheel brake* | 13 |
| 2.0.14. Accessories | 13 |

1.0. INTRODUCTION

DEAR CUSTOMER

Congratulations on your new VELA Medical Chair. This user manual contains useful information about positioning, operating and maintenance.

Important

Please read this user manual thoroughly and keep it for later use. This user manual must always accompany the chair. The user manual is also available on VELA Medical's website: www.vela-medical.com, where you also can find other relevant information about the chair.

If you have any questions, please contact VELA Medical's distributor in your country.

VELA Medical reserves the right to make changes to this manual.

VELA Medical

Gøteborgvej 8-12
DK-9200 Aalborg SV
Denmark
www.vela-medical.com

1.1. WARNINGS



| | |
|----------------|--|
| Use | In case of an allergic reaction, stop using the chair and consult a doctor |
| Use | Only use original spare parts from VELA Medical |
| Use | Only use the supplied charger for charging the chair |
| Use | Do not make any changes to the accessories that come with the chair |
| Use | Only use the chair for the intended purpose |
| Use | Always have the chair brake on when you sit down or get up |
| Use | Never use the chair on surfaces with level differences as the chair may become unstable |
| Use | Maintenance other than cleaning should always be performed by a professional |
| Use | Never service the chair while it is in use |
| Use | Always follow the recommended maintenance and cleaning instructions |
| Use | The chair's functions are intended for the user. Only the user or trained professionals should adjust the chair while the user is sitting in it |
| Use | Never use the electric height adjustment for more than two minutes at a time |
| Use | In the event of an electrical error or unintentional height adjustment, push both of the buttons for height adjustment simultaneously for 10 seconds |
| Use | Risk of shearing or squeezing fingers in moving parts. Therefore, never use the chair if someone is touching underneath the seat |
| Choking hazard | Choking hazard on small parts. Never allow children to play with the chair, as small parts are removable |
| Choking hazard | If the chair is equipped with electric height adjustment, cables will be present. Never let children play with these |
| Electric noise | Using cables other than the standard ones may increase electrical emission or de-crease electrical immunity |
| Electric noise | Avoid placing the chair close to appliances that emit electrical noise |
| Electric noise | The chair may only be used with accessories mentioned in this user manual, as elec-tro-magnetic interference may otherwise occur |
| Electric noise | If the chair moves unintentionally or does not work as expected, find another location for the chair |

1.2. TROUBLESHOOTING

| Mistake | Possible cause and solution |
|--|---|
| Electric height adjustment does not work | There is no power on the chair's battery. Charge the battery with the supplied charger |
| Manual height adjustment does not work | There is too much load or uneven load on the chair. Try to adjust the height without sitting in the chair. If this does not work, the gas cartridge is defective and should be replaced. Contact your VELA Medical supplier |
| The brake function does not work | The brake is not activated. See section 2.0.9. |
| The armrest seems loose | One of the hand screws might be loose. Tighten the hand screw |
| The chair rolls badly | Wheels are blocked. Clean the wheels for wire residue and the like. |

If the above does not solve your problem, you should contact your VELA Medical supplier.

1.3. VELA TANGO WITH ELECTRIC HEIGHT ADJUSTMENT

1.3.1. POWER INPUT

Voltage: 24Vdc

Current: 0.5A

Power: 12VA





1.3.2. EMISSION

A VELA Medical chair with electrical functions has been tested and complies with current requirements regarding electromagnetic radiation (EMC requirements). The chair is tested to the following specifications:

:: Emission class CISPR 11, Group 1, Class B

See more details at www.vela-medical.com

1.4. LABELS USED ON THE CHAIR

| | | |
|---------|--|--|
| Warning | Risk of shear and squeeze: Keep hands and fingers away from moving parts |  |
| Icon | Read "Instruction For Use" before use |  |
| Icon | Part used: Type B cf. IEC60601-1 covering the following areas of the chair: |  |
| Labels | 60 x 22,5 mm The last six digits are the product ID number: This label only applies to chairs with electrical functions: The chair is IP21 approved, which corresponds to protection against water droplets, as well as against touching dangerous parts with a finger. |  <p><Model name> <Model ref. no.> <Manufacturer> <Address> <Webpage></p> <p>CE MEDICO <Safe working load> <Production date/order no.></p> <p>UDI: (xx) xxxxxxxxxxxxxxxx (xx) xxxxxx</p> <p>Follow instructions for use</p> <p>EN60601-1 IP21</p> <p>Duty cycle: 10%, max. 2min/18min</p> <p>Blade fuse: 10A/32V</p> <p>Input: 24V== 0.5A</p> <p>Storage limits -20°C / 60°C</p> |

Labels on the chair should be read at a distance of approx. 40 cm under normal daylight conditions.

6 For labels related to the operation of the chair, please refer to section 2.0.

IN GENERAL

1.5. SAFETY

This chair bears the CE mark and approved under the model name "VELA Tango" and conforms to all applicable EU rules. The chair is manufactured by VELA Medical, Goeteborgvej 8-12, 9200 Aalborg SV, Denmark.

The VELA Tango 310E with electrical functions meets all safety requirements according to IEC 60601-1, including requirements for EMC according to IEC60601-1-2.

1.5.1. INTENDED USE

This device is intended to support a bariatric user in a sitting posture while performing indoor activities of daily living, including work activities.

- :: The chair is not built to stand on
- :: The wheel brake must always be activated when sitting down or standing up

All variants of VELA Tango 310 are tested and approved for a maximum load of 200 kg. For optimal use, the chair's recommended user weight is 130-160 kg.

All variants of VELA Tango 310E are tested and approved for a maximum load of 300 kg. For optimal use, the chair's recommended user weight is 130-300 kg.

1.5.2. TRANSPORT

When transporting the chair by car, airplane or any other means of transport, always use the chair brake. The chair is not approved for use as a passenger seat in any means of transportation and may not be used as such.

The chair is designed to be transported under the following environmental conditions:

Temperature: - 20°C to + 60°C

Humidity: 15-95%

Allow the chair to stand for 24 hours before use if it has been stored below - 10°C.

1.6. WARRANTY

VELA Medical grants a right of claim under applicable law in the country where the chair is acquired. The right of complaint applies only when using original spare parts and accessories, as well as customisations made by professionals.

Note: if an exclusive contractual agreement has been made directly with VELA Medical, this applies!

VELA Medical assumes no liability for damages to the product or the user caused by:

- :: Transportation
- :: Misuse
- :: Normal wear and tear
- :: Wrong use
- :: Use of spare parts and accessories that are not original VELA Medical products
- :: Adjustments made by unauthorised service technicians

1.7. UNPACKING THE CHAIR

The chair comes in environmentally friendly and recyclable packaging.

In case of defects or transport damage to the chair, please contact your VELA Medical supplier right away. In such cases, do not attempt to use the chair or to fix it yourself.

On delivery you have received:

- :: VELA Medical Chair
- :: Instructions for use
- :: Charger (only for chairs with electrical functions)

1.8. MAINTENANCE

We recommend an annual inspection of the chair to identify any defective or worn parts. All removable parts should be tightened regularly. You can find various guides on www.vela-medical.com.

1.8.1. WHEELS

We recommended you regularly check the wheels and remove any residual threads, etc.

1.8.2. FRAME

Clean with lukewarm water and if necessary, add some pH-neutral detergent.

1.8.3. PADDING

Fabric: Clean by chemical cleaning, with furniture cleaner or with water using a hard-wrung cloth.

Synthetic fabric or leather: Clean with water or mild soap using a hard-wrung cloth.

Medical textiles: Clean as described above. Can also be cleaned with disinfectants and alcohol. To see which specific concentrations the different pad types can handle, see www.vela-medical.com

Note: On www.vela-medical.com you can find cleaning instructions for all our standard upholstery types.

1.8.4. RECYCLING GUIDE

The chair's life expectancy is five years for regular and proper use. We recommend that VELA, your VELA supplier or a trained professional perform an annual inspection with the replacement of worn parts.

1.9. RECYCLING AND DISPOSAL

This product contains recyclable materials. Therefore, the product must be disposed of in accordance with local regulations and not together with ordinary household waste. Proper disposal and recycling will help prevent negative effects on the environment and people.

Note: Electrical and electronic equipment and accompanying batteries contain materials, components and substances that can be harmful to human health and the environment if the waste is not handled properly.


Electrical and electronic equipment and batteries are marked with the icon shown below. This symbolises that electrical and electronic equipment and batteries must not be disposed of together with unsorted household waste but must be collected separately.



It is important that you return your used batteries to the established schemes. In this way, you help ensure that batteries are recycled in accordance with current legislation and do not cause unnecessary harm to the environment.

2.0. USER MANUAL


The gas spring and the electric lifting mechanism should only be replaced by a VELA Medical technician or a person trained in this field.



2.0.1. HEIGHT ADJUSTMENT (EL)*


Can be operated by pressing on the switch that is located on the side of the right or left armrest.

Note: VELA Tango El is also available with height adjustment via handheld control (see item 2.0.3.).




2.0.2. HEIGHT ADJUSTMENT (MANUAL)*

Adjusted by pulling up the handle under the armrest on the right-hand side.




2.0.3. HANDHELD CONTROL (EL)*

The handheld control is used to adjust the height of the chair via push buttons on the front.




2.0.4. ADJUSTMENT OF SEAT DEPTH

The seat depth is adjusted by loosening the 2 Unbrako bolts under the seat (next to the backrest arm), after which you will be able to push the backrest back and forth. Remember to tighten the Unbrako bolts securely again after adjusting the seat depth.




2.0.5. ADJUSTING THE BACKREST ANGLE

The backrest is angled by loosening the small screw on the back bolt, after which the backrest can be adjusted to the desired angle. Remember to tighten the bolt securely once the correct angle has been achieved.




2.0.6. SETTING THE BACK HEIGHT

The back height is adjusted by loosening the 4 Unbrako bolts that are placed behind the backrest arm. Once they are loose, the height of the back can be adjusted. Remember to tighten the Unbrako bolts securely again.



2.0.7. ADJUSTING THE ARMREST HEIGHT

Adjust by loosening both screws under the armrest bar and adjusting the height of the armrest. Remember to tighten the screws securely again.




2.0.8. ADJUSTING THE ARMREST WIDTH

Adjust by loosening both screws under the seat. Remember to tighten the screws securely again.

Note: The armrests may only be adjusted up until the red 'stop' mark.


** Depending on model*



2.0.9. BRAKE*

The brake on the chair is activated by moving the handle on the side of the chair forward. Unlock the chair by moving the handle back into the upright position.


Note: Available with the brake mounted on the right or left side.



2.0.10. CHARGING (EL)*


Daily charging is recommended. The chair is charged by plugging the charger into the switch placed behind the chair under the seat on the right side. An LED will light yellow when the chair is charging and green when the chair is fully charged.

Note: The chair should not be used during charging.



2.0.11. STROLLING BRACKET*

Strolling bracket for transporting a person in the chair. The strolling bracket is fitted to the top back clamp.




2.0.12. ELECTRIC BRAKE (EL)*

Chairs with electric brake, is operated via remote control. We advice to fully brake/unlock the brake every time.

Note: This function can only be controlled via remote control.

* Depending on model



2.0.13. 4-WHEEL BRAKE*

Is activated via the hand/foot brake on all chairs except on chairs with electric brake, where it is activated via hand box or switches. 4-wheel brake brakes on all wheels and provides a stable seat on uneven/skewed surfaces.

Note: On the electric version it is advisable to activate or deactivate the brake completely when using this feature.

2.0.14. ACCESSORIES

This chair is available with a range of accessories. See a complete list of accessories at: www.vela.eu

Note: Charger is included as an accessory when purchasing a VELA chair with electrical functions.

* *Depending on model*

